MEDsan® SARS-CoV-2 Antigen Rapid Test

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Email: kommunikation@sanagroup.com

MEDsan® SARS-CoV-2 Antigen Rapid Test

MEDsan® SARS-CoV-2 Antigen Rapid Test is a solid phase immunochromatographic assay intended for the *in vitro* qualitative detection of specific severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) antigens in human nasopharyngeal and/or oropharyngeal secretion.

- High delivery capacity
- Highest quality standards
- High reproducibility

Advantages

- » Quick result - 15-20 minutes
- » Easy to use - everything from one set
- » Possibility to test whenever there is no laboratory available
- » Test for professional use only

Kit Equipment

- » Test cassette (in foil pouch including desiccant)
- » Extraction buffer
- » Extraction tube with dispensing tip
- » Specimen collection swab
- » Package insert (DE, EN, IT, SP, FR, PL, PT, RU)

Format:	Lateral flow test / Immunochromatographic <i>in vitro</i> test
Type:	Qualitative
Instruments:	No further devices required
Specimen type:	Nasopharyngeal or oropharyngeal swab
Target:	Antigen nucleocapsid (N)
Reading time:	15 – 20 minutes
Specificity:	99.8%
Sensitivity:	92.5%, for samples with C_t values $\leq 33 = 96.5\%$
Limit of detection (LoD):	14.4 TCID ₅₀ /mL
Storage temperature:	2 – 30 °C
Available at:	Pharmaceutical wholesaler, medical specialist trade, laboratory specialist trade
PZN, PPN:	16860851 / 111686085184 (Single box) 16860868 / 111686086874 (Box of 25)
Pack size:	Box of 25 / Single box including all components
Allowed users:	Healthcare professionals
Package insert:	8 languages (DE, EN, IT, SP, FR, PL, PT, RU)

Orderhotline: +49 40 38 02 38 38 | info@medsan.eu

MEDsan

biological health solutions

MEDsan® SARS-CoV-2 Antigen Rapid Test

Quick start guide


Figure 1:
Healthcare professional specimen collection

IVD CE


Figure 2:
Fill in extraction buffer


Figure 3:
Dissolve the sample in the extraction buffer by rotating


Figure 4:
Attach the dropper for precise dosing


Figure 5:
Apply 2 drops

15 minutes

Figure 6:
Read result within 15-20 minutes

Simple result control

A validated result requires the illustration of the control line "C".

Positive:


The result is positive as soon as the test line "T" appears.

Negative:

The result is negative if no band forms for test line "T".

Invalid:

The result without control line "C" is invalid.

Positive

Negative

Invalid

Performance:

MEDsan biological health solutions		RT PCR		
		Positive	Negative	Total
MEDsan® SARS-CoV-2 Antigen Rapid Test	Positive	111	1	112
	Negative	9	499	508
	Total	120	500	620
Sensitivity		92.5% [95% CI: 86.4%, 96.0%]		
Specificity		99.8% [95% CI: 98.9%, 100%]		
Accuracy		98.4% [95% CI: 97.1%, 99.1%]		

Limit of detection (LoD): 14.4 TCID₅₀/mL

More questions? Contact MEDsan®:

MEDsan GmbH

Brandshofer Deich 10 | 20539 Hamburg
Germany

Phone: +49 40 38 02 38 38

Fax: +49 40 228 58 56 99

Email: kommunikation@sanagroup.com